



## **PATHWAY STUDIO® DESKTOP 7 SYSTEM REQUIREMENTS**

### **OPERATING SYSTEM**

Windows XP Professional with Service Pack 2 (or later)  
Windows Vista (Home Basic, Home Premium, Ultimate, Business or Enterprise edition)  
Windows 7 (Home Basic, Home Premium, Ultimate, Professional or Enterprise edition)  
\*32-bit or 64-bit Edition supported, 64-bit Edition recommended.

### **RAM**

Windows XP - 2 GB Minimum (3 GB recommended)  
Windows Vista - 2 GB Minimum (4 GB or more is recommended)

### **MEDIA DRIVES**

DVD-ROM drive (for execution of Pathway Studio® installation media)

### **HARD DISK SPACE USAGE**

ResNet® 7 Mammalian database - 10 GB  
ResNet® Plant database - 3 GB  
Pathway Studio® application, including MedScan Reader - 1 GB  
MS SQL Server 2008, Workgroup Edition – 1 GB

### **OTHER REQUIREMENTS**

NTFS file system  
Microsoft Internet Explorer 6 or later

Microsoft .NET 3.5 SP1 (or 3.5.1) is required, and provided on the Pathway Studio® installation media or downloaded from Microsoft during installation.

Microsoft MS-SQL 2008 Workgroup Edition is required, and is provided on the Pathway Studio® installation media and automatically installed during Pathway Studio® installation. If the installation finds an ARIADNE instance of MS-SQL 2005 Standard Edition already installed, the ResNet 7 database will be installed into that, and installation of MS-SQL 2008 will not be initiated. Note: use of MS-SQL requires a license agreement with Microsoft which can be provided by Ariadne at a discount when licensing Pathway Studio®.