

PATHWAY STUDIO®

Software for Visualization and Analysis of Biological Pathways

PATHWAY STUDIO 5.0 DESKTOP SYSTEM REQUIREMENTS

OPERATING SYSTEM

Windows 2000 Professional Edition with Service Pack 4

Windows XP Professional with Service Pack 2

Windows Server 2003 (Enterprise, Standard or Datacenter edition) with Service Pack 1 or later

Windows Small Business Server 2003 with Service Pack 1 or later

Windows Vista (Ultimate, Business or Enterprise edition)

RAM

1 GB (2 GB is recommended)

MEDIA DRIVES

DVD-ROM drive

HARD DISK SPACE USAGE

| | |
|--|--------------|
| ResNet 5.0 database ¹ | 10 GB |
| Plant ResNet database | 3 GB |
| Pathway Studio Client with Medscan Reader module | 500 MB |
| MS SQL Server 2005 Standard edition | 1 GB |
| Overall HDD space | 15 GB |

OTHER REQUIREMENTS

Microsoft Internet Explorer 6.0 with Service Pack 2 or later

¹ The logical disk on which ResNet 5.0 is installed must be formatted with NTFS file system