

# PATHWAY STUDIO®

Software for Visualization and Analysis of Biological Pathways

## PATHWAY STUDIO DESKTOP 6 SYSTEM REQUIREMENTS

### OPERATING SYSTEM

Windows XP Professional with Service Pack 2 (or later)

Windows Vista (Home Basic, Home Premium, Ultimate, Business or Enterprise edition)

### RAM

Windows XP - 1 Gb (2 Gb is recommended)

Windows Vista - 2 Gb (3 Gb or more is recommended)

### MEDIA DRIVES

DVD-ROM drive

### HARD DISK SPACE USAGE

ResNet® Mammalian 6 database - 10 GB

ResNet Plants 2 database - 3 GB

Pathway Studio® application - 500 MB

MS SQL Server 2005, Standard Edition – 1 GB

### OTHER REQUIREMENTS

NTFS file system

Microsoft Internet Explorer 6.0 with Service Pack 2 (or later)